

AirSep® FreeStyle®

PORTABLE OXYGEN CONCENTRATOR

FREEDOM HAS A NEW LOOK

CAIRE introduces its original lightweight portable oxygen concentrator — the AirSep® FreeStyle®. The new design adds features and benefits that will not only make the device more appealing to patients, but even more reliable for providers.

- **New Power Supply Connection**

Completely re-designed connector featuring multi-directional installation with no pins to align and no release button to press is simpler than ever to use. Robust metal connection points ensure the highest durability of both the power supply and FreeStyle for patients and provider alike.

- **New User-Replaceable Battery**

Swap battery packs on the go instead of waiting to re-charge an internal battery or having to add heavy external batteries. The battery pack will re-charge anytime AC or DC power is connected and is easily accessible from the front of the device.

- **Longer Battery Life**

The new battery pack provides an additional 1 hour of duration at setting 2, providing up to 3 ½ hours with a single battery pack at setting 2.

- **New Color Scheme**

New color scheme is more appealing and non-medical in appearance. The matte finish will help resist everyday scratches.




CAIRE®

A Chart Industries Company

770.721.7759 | Toll Free 1.800.482.2473

www.CAIREMedical.com

FreeStyle Concentrator ON SALE



** SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.

PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, AND DETAILED SAFETY INFORMATION.

© 2016 CAIRE Inc. All Rights Reserved ML-CONC0078 Rev B

PROVIDER FEATURES

- Settings 1-3 cover the most common settings without the added weight or cost of larger capacity units
- Non-delivery modality that eliminates expensive cylinder deliveries.
- Time-proven and trusted quality
- Standard 3-year warranty includes sieve beds

PATIENT FEATURES

- Under 5 pounds with a small compact frame
- User-replaceable battery pack with up to 3½ hour battery life at setting 2
- Easily wearable with a handle, shoulder strap, or backpack
- FAA Approved
- Fully functional on AC or DC power
- Simple, One-Touch Operation

CLINICAL FEATURES

- UltraSense™ Technology provides sensitive breath detection
- Quickly detects and pulses oxygen to promote oxygen delivery in the earliest part of the breath
- Adjusts bolus size based on breath rate
- Lightweight package promotes greater ambulation and quality of life



SPECIFICATIONS

Flow Rates	Pulse 1-3
Dimensions	8.6 x 6.1 x 3.6 in (21.8 x 15.5 x 9.1 cm)
Weight (w/Battery)	4.9 lbs (2.2 kg)
AC Power	100-240 VAC, 50-60Hz (1.75 amp)
DC Power	11-16 VDC (8.0 amp max)
Operating Temp.	41°F to 104°F (5°C to 40°C)
Sound Level	41 dB(A) @ Setting 2
Operating Altitude	Up to 12,000 ft (3657.6 m)
O2 Concentration	90% (+5.5% / -3%)
Battery Duration	Setting 1 - 6.5 hours (6:30) Setting 2 - 3.5 hours (3:30) Setting 3 - 3.0 hours (3:00)
Battery Re-Charge	4 - 5 hours
Warm-Up Time	< 2 minutes
O2 Output	Setting 1 - 227.5 mL/min Setting 2 - 420 mL/min Setting 3 - 500 mL/min
Limited Warranty	3 years unit ** 1 year battery & power supplies

ORDERING INFORMATION

FreeStyle Includes: Portable Concentrator, Battery Pack, AC/DC Universal Power Supply, AC Power Cord, DC Power Cable, Carrying Case, Shoulder Strap, User Manual, Carry-all Accessory Bag

AS095-101	USA
AS095-102	European
AS095-B101	USA - 2 batteries
AS095-B102	Euro - 2 batteries
AS095-103	Australia
AS095-104	UK
AS095-B104	UK - 2 batteries
AS095-105	Brazil
AS095-106	Swiss

other countries' power supplies may be available - contact customer service

ACCESSORIES

BT033-1	Spare Battery Pack
MI284-1	Backpack Harness
MI419-1	Carrying Case
MI320-1	Carry-All Bag
MI426-1	Shoulder Strap
F1144-1	Intake Filter
PW032-1	Universal Power Supply (US)
PW032-2	Universal Power Supply (EU)
PW032-4	Universal Power Supply (UK)
<i>Power Supply Included</i>	
BT032-1	Desktop Charger (US)
BT032-2	Desktop Charger (EU)
BT032-4	Desktop Charger (UK)
<i>Cords Included</i>	


CAIRE®
A Chart Industries Company

770.721.7759 | Toll Free 1.800.482.2473

www.CAIREMedical.com



** SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.

PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, AND DETAILED SAFETY INFORMATION.

© 2016 CAIRE Inc. All Rights Reserved ML-CONC0078 Rev B