

**Affordable Price.  
Uncompromised Quality.**



[www.aedstore.in](http://www.aedstore.in)

# **AED Cabinet with Alarm**

**Model: CABBY-ECO**





www.aedstore.in

# Key Features

Model: CABBY-ECO



Easily Replaceable Button Battery



Alarm Mute Button



T-Handle Panel Lock



All three-side Red & Black AED Signage (Left, Right and Front)



Powder Coated Paint



4-Point Wall Mounting System



Polycarbonate Glass



[www.aedstore.in](http://www.aedstore.in)

# Universal Design Fit All AED's

Model: CABBY-ECO



Heartplus



BPL



Cardiac Science



Nihon Kohden



Schiller



Zoll



Philips



AED 7000



Radian



Defibtech



Physio-Control



Mediana





www.aedstore.in

# Specifications

**Model: CABBY-ECO**

<b>Physical</b>	
<b>Dimension</b>	(W) 39.2cm x (H) 42cm x (D) 17cm
<b>Weight</b>	5.200 Kg +/- 100gms
<b>Construction</b>	Made from Heavy duty CRC Steel Sheet
<b>Paint</b>	Oven baked Powder coated paint
<b>Paint Color</b>	Red, White
<b>Locking</b>	Safety Lock with Two keys
<b>Glass Window</b>	Clear PC glass window with red & black lettering to facilitate visual identification of the AED
<b>AED Signage</b>	All three-side red & black AED Signage (Left, Right and Front)
<b>Wall Mounting</b>	Four holes are punched in the back of this cabinet for a wall mount

<b>Battery</b>	
<b>Type</b>	Button Battery
<b>Quantity</b>	3 Cell
<b>Standby Life</b>	More than 1 year

<b>Fit to Brands</b>	
<b>Brand of AED</b>	All Brands Schiller, Nihon Koden, Cardiac Science, Schiller, Zoll, BPL, Defibtech, Physio-Control, Radian, HeartPlus, Philips, Mediana, M&B etc.
<b>Customization</b>	Yes Possible
<b>Branding</b>	ODM is possible

<b>Features &amp; Safety</b>	
<b>Alarm</b>	Cabinet is enables with 80 dB Alarm
<b>Sound</b>	80 dB
<b>Protection</b>	Protection of AED from Dirt, Dust and possible theft
<b>Warranty</b>	Two years



+91 9871181879

## Essential Medical

Address : B-14, Durga industrial Park Prakash Industrial Estate,  
Sahibabad Industrial Area Site 4, Ghaziabad, U.P.- 201005  
E-mail : [ajay@akscorp.org](mailto:ajay@akscorp.org)  
Website: [www.aedstore.in](http://www.aedstore.in)