



F4 Full Face Mask

ULTRA-LIGHTWEIGHT AND SUPERIOR COMFORT

The full face mask single-layer soft cushion is specially designed to eliminate bruising and create a perfect fit to the contours of the face. The Full Face Mask (F4) channels airflow non-invasively to a patient from a positive airway pressure device such as a continuous positive airway pressure (CPAP) or bi-level system (BiPAP).



AutoFit cushion

The shaped groove on the cushion provides an inclusive range adapting to most of the face contours.

Slim forehead brace

A balance between the reinforced stability and the widened vision is achieved by the slim brace inspired by the minimalist concept.



Multifunctional port

Two embedded ports on the bilateral sides of the main frame make the oxygen input or pressure test as convenient as possible.

Ball & socket elbow

The universal elbow allows the user to turn over freely as needed.



BMC Medical Co. Ltd.

Room 110, Tower A Fengyu Building, No.115 Fucheng Road,
Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF CHINA

Tel: +86-10-5166 3880

Fax: +86-10-5166 3880 ext. 810

en.bmc-medical.com

SC-02-21-06-004-Ver2.1

