



EVOX

Contributing To Healthy Life

CPAP / BIPAP SERIES



For Sleep Apnea Hypopnea
Syndrome(SAHS)



Technical Data



	CPAP-C09	BIPAP-B19
	C09	B19
Mode		
CPAP	✓	✓
Auto CPAP	✓	
S		✓
T		✓
S/T		✓
Parameters		
Pressure Range	4-20 cmH2O	4-25cmH2O
Belex	2-4cmH2O	2-4cmH2O
Ramp	0-60 min	0-60 min
Esens (1-6 levels)		✓
Isens (1-6levels)		✓
Auto-leak Compensation	✓	✓
Auto-altitude Adjustment	✓	✓
Auto-On/Off	✓	✓
Pressure Titration	✓	
Noise Level (10cmH2O)	<30dBA	<35dBA
Heating Humidifier		
Adjustable Humidify	✓	✓
Replaceable Water Chamber	✓	✓
Data Download		
Data Details	✓	✓
Data Summary	✓	✓
Download to PC	✓	✓