

Rad-G™ Pulse Oximeter


A rugged handheld device powered by clinically proven Masimo SET® pulse oximetry



Strong external casing built to withstand drops from 6 ft. onto a hard surface



Light weight and slim profile allow for easy transport and storage



Internal rechargeable battery provides 24 hours of operational use



Reliable Monitoring

Hospital-grade Masimo SET® Measure-through Motion and Low Perfusion™ pulse oximetry*

Informative Display

Display of pleth waveform for additional information to use in patient assessment

Alarm Management

Configurable alarms to offer insight into patient status changes that may require intervention



Available with Masimo single-patient-use and reusable sensors, Rad-G monitors:



Clinically Proven Technology Used in Leading Hospitals Worldwide

per year around the world¹

hospitals according to the 2020-2021
U.S. News & World Report
Best Hospitals Honor Roll²

have shown that Masimo SET[®]
outperforms other pulse
oximetry technologies³