

BIPAP 25 & 30 PLUS




MODELS	BPAP25PLUS	BPAP30 PLUS
Modes	S, T, S/T and CPAP	S, T, S/T and CPAP
IPAP Range	4-25 CmH2O	4-30 CmH2O
EPAP range	4-20 CmH2O	4-20 CmH2O
CPAP range	4-20 CmH2O	4-20 CmH2O
1-trigger	Auto,1-6	Auto,1-6
E-trigger	Auto,1-6	Auto,1-6
IPAP rise time	Auto,0.1-0.6s	Auto,0.1-0.6s
Breath rate	2-40 bpm	2-40 bpm
Noise Level (dB)	<28dB	<28dB
Ramp time(minute)	0-60	0-60
Display data	IPAP, EPAP, Rate(bpm), LE, Vt. Vleak, Mvent, Trigger 1/E	IPAP, EPAP, Rate(bpm), LE, Vt. Vleak, Mvent, Trigger 1/E
Alarms	Tube Alarm, Apnea Alarm, High/low rate Alarms & power failure	Tube Alarm, Apnea Alarm, High/low rate Alarms Mvent alarm & power failure
Dimensions	18*27*11cm	18*27*11cm
Weight	1.7 Kg	1.7 Kg
Screen size	4.3 inch	4.3 inch
Heated humidifier	Integrated	Integrated
Memory	Micro SD Card 8 GB	Micro SD Card 8 GB
Warranty	2 years	2 years

BIPAP 25 & 30 PLUS



Electrical Characteristics	
Electrical supply	100-240V, 50-60 Hz
Power Consumption	<80W with Humidifier
Environmental requirements	
Temperature	5°C-35°C
Humidity	<80W with Humidifier
93/42/EC Classification	30-95%, no condensation
Altitude	From Sea level to 2700 mtrs
Storage	-20-60-°C



Mode->S/T	
Mode	S/T
IPAP	4.0 cmH ₂ O
EPAP	0.40 cmH ₂ O
Rise Time	0.1 sec
Rate	10.0 bpm
Insp Time	1.5 sec
I-trigger	6